



### BONE BRIDGE GUIDE

<b>Warnings / Precautions</b>	Before use, carefully inspect the burr for any irregularities or eccentricities. Damaged or worn product should be replaced Cold soak in glutaraldehyde, chlorine, or ammonium solutions, or dry heat sterilization is not recommended as damage to the instrument finish may occur.			
<b>Limitations</b>	Discard any device that shows signs of damage or wear			
INSTRUCTIONS				
<b>Point of Use</b>	This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse. Promptly and thoroughly rinse with deionized water after each use.			
<b>Containment and Transportation</b>	It is recommended that devices are reprocessed as soon as is practical following use.			
<b>Preparation for decontamination</b>	Promptly and thoroughly rinse with deionized water after each use.			
<b>Cleaning: Automated (Do NOT use ultrasonic washer)</b>	Remove instruments and equipment from any sterilization trays before placing into washer baskets. Orient devices following recommendations of washer/disinfector manufacturers. Use alkaline or neutral pH detergent recommended by washer/disinfector or detergent manufacturers. These products have been validated for effective cleaning using an automatic washer/disinfector cycle consisting of a minimum 44 minutes total time, including a pre-wash, main wash & rinse, and thermal rinse. The thermal rinse shall be at least 10 minutes long at a minimum temperature of 60°C.			
<b>Cleaning: Manual</b>	Soak in lukewarm*, mild* enzymatic detergent, and deionized water for a minimum of two minutes. Then clean ultrasonically in lukewarm* solution of mild* detergent and deionized water for at least 30 seconds. Rinse thoroughly with deionized water and wipe dry. *less than 43°C; pH 7.0 - 8.5 Note: When using an ultrasonic cleaner or a spray washing machine, follow the manufacturer's recommendations, particularly with regard to articulated instruments and positioning of instruments.			
<b>Disinfection</b>	No particular requirements			
<b>Packaging</b>	No particular requirements			
<b>Sterilization (Temperatures are minimum required, times are minimum required)</b>	<b>Cycle:</b>	<b>Gravity</b>	<b>Pre-vac</b>	
	<b>Temperature:</b>	121°C	134°C	
	<b>Time:</b>	30 min 8 min vacuum	3 min 8 min vacuum	
	<b>Drying: 8 min or until visibly dry</b>			
	<b>100% EtO</b> Not validated		<b>STERRAD Sterilization</b> Not validated	
<b>Maintenance, Inspection and Testing</b>	No particular requirements. Discard any device that shows signs of damage or wear			
<b>Storage</b>	No particular requirements			
<b>Additional Information</b>	None			

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieve the desired result. This normally requires validation and routine monitoring of the process.