



ENDOSCRUB® 2 FINGER SWITCH

Warnings / Precautions	Disconnect the finger switch from the EndoScrub 2 pump before cleaning.			
Limitations	After cleaning and sterilization, verify functionality prior to re-use.			
INSTRUCTIONS				
Point of Use	This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse. After use, remove the finger switch from the sheath and disconnect the plug from the pump. Thoroughly rinse with water following use.			
Containment and Transportation	No particular requirements. It is recommended that devices are reprocessed as soon as is practical following use.			
Preparation for decontamination	No particular requirements.			
Cleaning: Automated	Not validated			
Cleaning: Manual	Dip the finger switch housing in a diluted mixture of mild (pH 7.0 - 8.5) enzymatic detergent. (Follow detergent manufacturer's instructions for proper dilution.) Thoroughly clean the housing with a soft instrument brush to remove any blood and tissue. Rinse the housing thoroughly with tap water and wipe dry. Note: If wiping the cord dry, be sure to hold the cord and not the housing to avoid stressing or breaking the electrical connections located inside the housing.			
Disinfection	No particular requirements.			
Packaging	No particular requirements.			
Sterilization (Temperatures are minimum required, times are minimum required)	Cycle:	Gravity	Pre-vac	
	Temperature:	132 °C	132 °C	
	Time:	10 min	10 min	
	Drying: 8 minutes, or until visibly dry			
	100% EtO Not validated		STERRAD Sterilization Not validated	
Maintenance, Inspection and Testing	Inspect finger switch for any damage before and after each use. If damage is observed do not use the finger switch until it is repaired or replaced. After cleaning and sterilization, verify functionality prior to re-use.			
Storage	No particular requirements.			
Additional Information	None.			

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieve the desired result. This normally requires validation and routine monitoring of the process.