



ORO-TRACHEAL INJECTION INSTRUMENT

Warnings / Precautions	Hypodermic needles used with this injection instrument are intended for single use only, do not resterilize or reuse. Cold soak in glutaraldehyde, chlorine, or ammonium solutions, or dry heat sterilization is not recommended as damage to the block may occur.				
Limitations	After cleaning and sterilization, verify functionality prior to re-use.				
INSTRUCTIONS					
Point of Use	This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse. Rinse immediately after use in distilled or demineralized water				
Containment and Transportation	No particular requirements It is recommended that instruments are reprocessed as soon as is practical following use.				
Preparation for decontamination	No particular requirements				
Cleaning: Automated (Do NOT use ultrasonic washer)	Remove instruments and equipment from any sterilization trays before placing into washer baskets. Orient devices following recommendations of washer/disinfector manufacturers. Use alkaline or neutral pH detergent recommended by washer/disinfector or detergent manufacturers. These products have been validated for effective cleaning using an automatic washer/disinfector cycle consisting of a minimum 44 minutes total time, including a pre-wash, main wash & rinse, and thermal rinse. The thermal rinse shall be at least 10 minutes long at a minimum temperature of 60°C. Instruments with moveable parts should be lightly lubricated with a water soluble instrument (i.e. instrument milk) after each cleaning and before sterilization, to assure proper function. Do not use Silicone Spray.				
Cleaning: Manual	Disassemble instruments as necessary. Soak in lukewarm (less than 43°C), mild (pH 7.0 - 8.5), enzymatic detergent and deionized water for a minimum of 2 minutes, then clean ultrasonically in a lukewarm (less than 43°C), mild (pH 7.0 - 8.5), detergent and deionized water for at least 30 seconds. Instruments with moveable parts should be lightly lubricated with a water soluble instrument (i.e. instrument milk) after each cleaning and before sterilization, to assure proper function. Do not use Silicone Spray. Note: When using an ultrasonic cleaner or a spray washing machine, follow the manufacturer's recommendations, particularly with regard to articulated instruments and positioning of instruments.				
Disinfection	No particular requirements				
Packaging	No particular requirements				
Sterilization (Temperatures are minimum required, times are minimum required)	Cycle:	Gravity			
	Temperature:	121°C			
	Time:	30 min			
	Drying:	8 min or until visibly dry			
	EtO: Not validated		STERRAD Sterilization: Not compatible		
Maintenance, Inspection and Testing	After cleaning and sterilization, verify functionality prior to re-use.				
Storage	No particular requirements				
Additional Information	None				

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieve the desired result. This normally requires validation and routine monitoring of the process.